

# ACQ-06 Column Switcher

## —Provide Fast Purification

ACQ-06 column switcher allows chemists to switch automatically from one column to another, without evaporating solvents. It also offers multi-stationary phases for a complete chromatography solution.

- Time Saving  
Start a new run by switching to a new column
- Efficient Purification  
Further purification can be achieved through tandemed columns
- Good Compatibility  
Compatible with the purification system from most of the manufacturers in the market




# Structure



- ① Indicator Lamps  
Observe the status of the channels
- ② Column Support  
Suitable for columns of different specifications  
Enable column tandem
- ③ Standard Adaptor  
Compatible with flash columns from other manufacturers
- ④ Channel Switcher  
Avoid cross contamination

# Software





- Easy to Operation  
Six channels enable the switching from column 1 to column 6
- Independent Software  
Allow single use or connected with CHEETAH program

# Still purify your sample by hand?

Free yourself by choosing a fully automatic purification kit !

Just build a connection between

- ATS Auto-Sampler
- ACQ-06 Column Switcher
- CHEETAH™ System



MZ-Analysentechnik GmbH  
Wöhlerstraße 2-6 • D-55120 Mainz  
Tel +49 6131 68 66 19  
Fax +49 6131 68 66 20  
e-mail: info@mz-at.de  
www.mz-at.de

ABN: ZL-03210  
2012.10. Rev.1  
©Bonna-Agela Technologies  
All rights reserved.