


Chiral Technologies

USP List of HPLC Columns

USP Code Column Name	USP Specification Analysis (Pending/Official)	Drug Chemical Name	Drug Brand Name
L32 CHIRALPAK® WH	Limit of Organotin	Fluorodopa	FDOPA
L40 CHIRALCEL® OD CHIRALCEL OD-H	Limit of efavirenz R-enantiomer Limit of voriconazole-related compound B Limit of atomoxetine S- enantiomer Limit of valsartan-related compound A Limit of irinotecan hydrochloride enantiom	Efavirenz Voriconazole Atomoxetine Hydrochloride Valsartan Irinotecan Hydrochloride	Sustiva VFEND Strattera Diovan Camptosar
L41 CHIRALPAK AGP	Limit of ropivacaine-related compound B Limit of galantamine 4R, 8R-stereoisomer Limit of montelukast S-enantiomer Limit of tenofovir disoproxil-related compo	Ropivacaine Galantamine Hydrobromide Montelukast Sodium Tenofovir Disoproxil Fumarate	Naropin Razadyne Singulair Viread
L51 CHIRALPAK AD CHIRALPAK AD-H CHIRALPAK AD-3	Limit of abacavir enantiomer Limit of levetiracetam R-enantiomer Limit of tolterodine S-enantiomer Limit of atorvastatin 3S,5S-steriomer Measurement of diastereoisomeric ratio Limit of paroxetine-related compound C Limit of brinzolamide-related compound A	Abacavir Sulfate Levetiracetam Tolterodine Tartrate Atorvastatin Calcium Fulvestrant Paroxetine Hydrochloride Brinzolamide	Ziagen Keppra Detrol Lipitor Faslodex Paxil Azopt
L66 CROWNPAK® CR(+)	Limit of valacyclovir D-isomer	Valacyclovir Hydrochloride	Valtrex
L70 CHIRALCEL OC-H	Limit of oxaliplatin-related compound D	Oxaliplatin	Eloxatin
L79 CHIRALPAK HSA	Enantiomeric Purity	Calcium L-5- Methyltetrahydrofolate	Calcium Supplements
L80 CHIRALCEL OJ CHIRALCEL OJ-H CHIRALCEL OJ-3	Limit of Exemestane To be specified	Exemestane Clopidogral Bisulfate	Aromasin Plavix
L90 CHIRALPAK AS CHIRALPAK AS-H CHIRALPAK AS-3	Enantiomeric Purity	Frovatriptan Succinate	Frova

USP Code Column Name	USP Specification Analysis (Pending/Official)	Drug Chemical Name	Drug Brand Name
L93 CHIRALCEL OD-RH CHIRALCEL OD-3R	To Be Specified	Ezetimibe	Zetia
L99 CHIRALPAK IA CHIRALPAK IA-3	To Be Specified	Darifenacin Hydrobromide	Enablex
L107 CHIRALCEL OJ-RH	To Be Specified	Rosuvastatin Calcium	Crestor
L## (TO BE DECIDED) CHIRALPAK IC CHIRALPAC IC-3	Limit of atomoxetine S- enantiomer	Atomoxetine Hydrochloride	Strattera
L## (TO BE DECIDED) CHIRALCEL OD-R	To Be Specified	Latanoprost	Xalatan



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Wöhlerstraße 2-6 • D-55120 Mainz
 Tel +49 6131 68 66 19, Fax +49 6131 68 66 20
 e-mail: info@mz-at.de, www.mz-at.de

800 North Five Points Road, West Chester, PA 19380

800-6CHIRAL 610 594 2100 610 594 2325 FAX chiral@chiraltech.com www.chiraltech.com

