

Intrada Amino Acid

50 x 3 mm

Application

55 Amino acids (standard)

55種アミノ酸(標準試料)

Intrada Amino Acid, 50 x 3 mm

A: acetonitrile / tetrahydrofuran / 25mM ammonium formate / formic acid = 9 / 75 / 16 / 0.3

B: acetonitrile / 100mM ammonium formate = 20 / 80

0 %B (0-2.5min)

0-17 %B (2.5-6.5 min)

100 %B (6.5-10 min)

0.6 mL/min (6MPa), 35 deg.C, 5 uL (1000nmol/mL)

ESI (SIM, positive)

