



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
Wöhlerstraße 2-6 • D-55120 Mainz
Tel +49 6131 68 66 19
Fax +49 6131 68 66 20
e-mail: info@mz-at.de
www.mz-at.de

Transgenomic RPSep Fluoro Silica Columns (pentafluorophenylpropyl phase)

Finally! An improved Silica-based line of Reverse Phase Columns for those who need increased selectivity and reproducibility in their assays. Transgenomic is introducing special reverse phase column types to enhance your everyday separation needs, as well as to provide more power towards R&D projects.

Let us help you choose the perfect column for your separation needs

Benefits

- LC-MS Compatible
- Can be used with Reversed or Normal Phase and HILIC Mode
- Special Selectivity compared to Normal C18 phase

Uses

- Pharmaceuticals, natural products, research & method development
- Strong interaction with Basic, Halogenated analytes and Aromatics such as Fluorines and Amines



Beta Blockers, Fluorinated Corticosteroids, Nucleosides, Nucleotides, Purines, Pyrimidines, Taxanes and Tricyclic Antidepressants

Cat No.	Sphere	Length	ID		Cat No.	Sphere	Length	ID
SIL-99-1101	3µm	50mm	2.1mm		SIL-99-1106	5µm	50mm	2.1mm
SIL-99-1102	3µm	100mm	2.1mm		SIL-99-1107	5µm	100mm	2.1mm
SIL-99-1103	3µm	50mm	4.6mm		SIL-99-1108	5µm	50mm	4.6mm
SIL-99-1104	3µm	100mm	4.6mm		SIL-99-1109	5µm	100mm	4.6mm

Custom column dimensions also available



Corporate Headquarters:
Transgenomic, Inc.
12325 Emmet Street, Omaha, NE 68164
phone: +1(402) 452.5400 Fax: (402) 452 5401
email: info@transgenomic.com
On the web at: www.transgenomic.com

NEW

Powerful